

TERUFUSION BLOOD ADMINISTRATION SET

PRODUCT features	SPECIFICATION	BENEFITS
Spike	54.2+ /-5 (mm) length	Allows easy insertion into Transfusion port.
Protector	47.5+/-5 (mm) length	Tightly fitted in the spike. So sterility is ensured even after taking out from the pouch.
Big Double Chamber	140+ /-5 (mm) length	Allows smooth filtration and flow of blood.
Roller Clamp	45+/- 2 (mm) length	Gives perfect grip and flow can be controlled properly.
Filter Mesh	83 mesh (60mm length)	Ensures smooth filtration. No block during filtration.
Material used in Drip Chamber	Polypropylene	Bio-compatible
Sterility assurance level	(SAL) 10 ⁶	Ensured sterility
Quality system	EN ISO13485 : 2003	Assured as per standards
EOG Sterilization/ Validation & daily control	EN ISO 11135 - 1:2007	Assured as per standards
Biological Tests done (Sterility test done)	Pyrogenic Test Acute Toxicity Test Intradermal Reaction Skin Sensitivity Test Cell Toxicity Test Endotoxin Test	Assured sterility